

REPROKON MILK PROGESTERONE ENZYMEIMMUNOASSAY KIT



- **PROGESTERONE DETECTION**
- **EARLY PREGNANCY DIAGNOSIS**
- **OESRUS DETECTION**
- **TIMING OF INSEMINATION**
- **SILENT HEAT DETECTION**

TAEA-ANTHAM Progesterone EIA kit, "REPROKON" is used for the Progesterone measurements in total milk in dairy cows and based upon a simple, reliable and sensitive EIA technique. The range of total milk Progesterone concentrations used in the tests is 0-20 ng/ml. One kit contains six standards and enough chemicals for 40 duplicate samples, however, 80 samples could be analyzed by experienced persons by using them as singles.

QUALITY ASSURED N-15 ANALYSIS OF ENVIRONMENTAL SAMPLES

N-15 Laboratory in the Ankara Nuclear Agriculture and Animal Research Center is analyzing total nitrogen content and N-15 atom excess of domestic and foreign plant, soil and water samples as certified Laboratory under the Proficiency Testing Exercise, EQA2004, organized by IAEA Seibersdorf Laboratory.

For Information :
Turkish Atomic Energy Authority
Ankara Nuclear Research Center
in Agriculture and Animal Sciences
Tel: +90 312 815 43 00
Fax: +90 312 815 43 07

